

Plan Beyond.

Interventional Systems
B | BRAUN
SHARING EXPERTISE

ONOCOR
Advancing Safety First


ONOCOR[®] Retrieval Device
Transcatheter Foreign Object Retrieval

ŌNŌ® Retrieval Device: Key Attributes

The ŌNŌ Retrieval Device is designed to receive, align, compress and remove material from the vascular system, allowing clinicians to retrieve what previously may not have been in reach.^{1;2}

Designed to utilize the existing access site

Indicated for use in the cardiovascular system to retrieve foreign objects using minimally invasive surgical procedures²



The outside pull force translates into radial compression of the foreign object by the ŌNŌ basket, reducing the foreign object back to (or close to) the delivery profile and allowing for easy extraction³

In recent decades, what is achievable by the Interventionalist has advanced in complexity and scope.

The ŌNŌ Retrieval Device is an enabling technology designed to allow clinicians to plan beyond what was previously in reach. The ŌNŌ device is designed to receive, align, compress and remove material from the vascular system. It is compatible with commercially available vascular sheaths, endovascular snares and other graspers.

On the creative medical frontier, the ŌNŌ device is a must-have tool, complementing your skillset.

Receive. Align. Compress. Remove.



Tear-away sheath allows introduction into any retrieval sheath 12 F or greater (maximum 75 cm overall length)


Flexible trackable catheter with 7 F compatible lumen or smaller (minimum 100 cm length)

Proprietary 35 mm basket flare* compatible with vessels as small as 10 mm

Facilitates coaxial alignment of the object for streamlined capture into the retrieval sheath¹

Designed to enhance snaring capability and allows large device retrieval through a smaller sheath compared to snaring alone¹

Designed to allow unobstructed blood flow when deployed



* Note: Diameter of the basket barrel is 30 mm. Full diameter including the basket flare is 35 mm. Snare is not included.

Introduction and Burden

Intravascular or intracardiac foreign bodies, although rare, can cause severe morbidity.⁴



The continued advancement of medical technology

has led to the expansion of intravascular and intracardiac devices to treat a wide range of cardiovascular conditions.^{5; 6}



These include structural heart devices

(atrial septal defect [ASD] closure, left atrial appendage closure [LAA] devices, valves etc.), leadless pacemakers, coils, clips, balloons, guidewires, catheters, filters, and more.^{4; 5; 7}



As the use of these devices has increased,

so have the reports of complications such as device migration or embolization.^{8; 9}



Although some patients may be asymptomatic,

clinical signs can include arrhythmia, pulmonary symptoms, cardiogenic shock, and sudden cardiac arrest.^{7; 10}



Embolized devices or endovascular foreign bodies

often require urgent removal, using surgical or percutaneous techniques.^{4; 5; 7}



Device migration or embolization

can be immediate, taking place during the implantation procedure, or delayed, occurring within the hours, weeks, or sometimes months following the procedure.^{7; 11}



While uncommon, these events can be potentially life-threatening,

particularly if localized to the heart or great vessels.⁷



Balloons are often used in endovascular procedures

and can rupture which may lead to fragmentation and embolization.^{24; 25}

The incidence of balloon rupture is estimated to be between 3.6% to 10% for balloon angioplasty,²⁴ and up to 33% for balloon aortic valvuloplasty.²⁵ Balloon rupture does not necessarily lead to fragment embolization.

Reports have generally estimated the incidence of intravascular foreign body mispositioning or embolization to be around 0.1%–0.9%.⁴

However, rates may vary depending on the specific device, patient population, and procedure being conducted.

Device	Reported Device Embolization Migration or Dislodgement Rates*
Transcatheter edge-to-edge repair (TEER)/Clips	0.1% ^{12; 13}
ASD Closure Devices	0.1%–1.3% ^{7; 14; 15}
LAA Closure Devices	0%–2% ^{7; 16; 17}
Venous Stents	0.17%–4.3% ^{18; 19}
Leadless Pacemakers	0.13%–0.76% ^{20–23}

* Rates may vary depending on the specific device, patient population, and procedure.



Cardiac tumors and intravascular thrombi may migrate or localize in the heart requiring retrieval and removal.²⁶⁻²⁸


Intracardiac and intravascular thrombi are associated with significant morbidity and mortality.^{29; 30}

- Widely used management strategies include anticoagulation therapy, systemic thrombolysis, and surgical embolectomy.³¹
- Thrombi that are older, more organized and fibrin-rich often exhibit increased resistance to thrombolysis compared to acute red blood cell-rich thrombi.^{32; 33}
- In patients with catheter-related right atrial thrombi, primary anticoagulation and thrombolysis treatment fails in approximately 27.3% and 33.3% of cases, respectively.³⁴

Cardiac tumors or masses occur rarely.³⁵

- Primary cardiac tumors have an incidence of **0.002% to 0.3%** based on an autopsy series.³⁶
- The incidence of patients clinically diagnosed with a primary cardiac tumors, however, is estimated at 1,380 per 100 million individuals, suggesting that **many go undiagnosed.**³⁵
- Secondary cardiac tumors or cardiac metastases are approximately **22 to 132 times more common** than primary cardiac tumors.³⁵
- Surgical resection of cardiac tumors is traditionally conducted via open median sternotomy.³⁷

Recent reports have suggested that some cardiac tumors and intravascular thrombi can be successfully removed using less invasive or endovascular methods.^{27-29; 37}



Intracardiac devices such as leadless pacemakers may undergo planned retrieval and removal due to device malfunction or end of life.^{38; 39}

- The battery life of a leadless pacemaker is between 4.7 to 15 years depending on pacing thresholds. Given this, some patients may require the implantation of multiple devices throughout their lifetime.³⁹
- Simply abandoning the previous leadless pacemaker and implanting a new device may lead to compromised cardiac function, device-device interference, or other potential long-term adverse events.^{38; 39}
- Retrieving and removing the leadless pacemaker when possible, prior to implanting a new device would avoid these potential complications and reduce the amount of nonfunctional cardiac hardware present.^{38; 39}
- Successful and safe percutaneous removal of leadless pacemakers **4 to 5 years** after implantation have been reported, indicating that removal may remain feasible at later time points, such as at device end of life.⁴⁰



Clinical Value

The ŌNŌ Retrieval Device provides an effective percutaneous solution to retrieve objects from the cardiovascular system, avoiding the need for invasive open surgical removal procedures.⁴¹



The ŌNŌ Retrieval Device facilitates

percutaneous cardiac tumor removal before symptoms or complications develop, in patients where surgical resection is not yet indicated.³



Failure to remove intravascular foreign bodies percutaneously

may require revision to emergent surgical removal.^{4; 7; 42}



Surgical structural heart procedures

are associated with higher in-hospital mortality, greater complication risk, and longer length of hospital stay compared to transcatheter procedures.^{43; 44}



Emergency open surgical procedures due to transcatheter device complications

exhibit greater mortality than primary surgical structural heart procedures alone.⁴⁵⁻⁴⁷

The ŌNŌ Retrieval Device can be used to facilitate the capture and removal of a wide variety of intravascular foreign bodies with a high rate of clinical success, exhibiting considerable versatility and reliability.^{41; 48}

Intracardiac and intravascular thrombi are associated with significant morbidity and mortality.^{29; 30}






- The ŌNŌ Retrieval Device has been successfully used to retrieve atrial tumors, cardiac thrombi, ruptured balloons and associated fragments, leadless pacemakers, stents, clips, ASD devices, LAA devices, and more.^{11; 26-28; 48-51}
- A wide variety of other devices or objects may also be suitable for retrieval with the ŌNŌ Retrieval Device, including but not limited to ventricular septal defect (VSD) devices, patent ductus arteriosus (PDA) devices and other malpositioned or embolized occlusion devices, guidewires, catheter fragments, coils and more.⁴¹
- In addition to migrated or embolized foreign bodies, the ŌNŌ Retrieval Device may also be used to facilitate the removal of malfunctioning or non-functioning devices, including those at end of life such as leadless pacemakers.^{39; 41; 52}
- The ŌNŌ Retrieval Device is suitable for use during both emergency and planned intravascular retrieval procedures.⁴
- The ŌNŌ Retrieval Device has been **successfully deployed in 100% of clinical cases** to date.⁴⁸
 - Unsuccessful retrieval procedures were complicated by factors such as catheter sheath length, snare or grasper capture, and/or patient characteristics.¹



The ŌNŌ® Retrieval Device simplifies complex intravascular foreign body retrieval procedures and can easily be integrated into an interventional cardiologists' practice.^{3; 41}

- The ŌNŌ Retrieval Device assists with proper orientation or alignment of the object for appropriate capture and removal, such as with leadless pacemakers and clips where coaxial orientation within the sheath is required for removal.^{11; 51}
- The ŌNŌ Retrieval Device compresses the captured item, facilitating the potential removal of larger foreign bodies, or the use of smaller catheter sheaths.^{27; 41}
- Once the item is captured, the ŌNŌ Retrieval Device covers or shields it from critical surrounding cardiac or vascular structures, such as cardiac valves, protecting against potential trauma or damage.^{11; 50}
- The nitinol ŌNŌ basket is easily collapsed and redeployed during the procedure to optimize basket size, engagement angle, device alignment, or to withdraw the object to a safer location where appropriate alignment and compression for removal can take place without risk to critical structures.^{11; 50}
- The ŌNŌ Retrieval Device is compatible with a variety of commercially available vascular sheaths, endovascular snares, and other graspers, with no other specific accompanying equipment required, allowing clinicians continued use of their preferred retrieval tools.^{1; 3}
- Minimal training is likely required to use the ŌNŌ Retrieval Device, especially for clinicians familiar with catheter-based therapies.³ Continued expert teleconferencing or in-person support on use of the device is available to clinicians as needed.





The ONO Retrieval Device coaxially aligns the object for easy extraction from the body. Images provided by Onocor LLC.

Resource Use and Economic Value of the ŌNŌ® Retrieval Device

Percutaneous foreign body removal procedures using the ŌNŌ Retrieval Device may shorten procedure time compared to open surgical procedures.



A case report

of percutaneous resection of an atrial tumor using electrosurgery and the ŌNŌ Retrieval Device had a recorded procedure time of **43 minutes**.²⁸



The average operation time

for surgical cardiac tumor removal is **198 minutes**.⁵³



Endovascular foreign body removal procedure time

is estimated to range from **6 to 208 minutes**.^{5; 54}



The mean procedure time

of mitral valve surgery via conventional midline sternotomy is **210.7 minutes**.⁵⁵




TAVR and PCI procedures

require on average less time compared to SAVR and CABG procedures (**54.6 minutes vs. 239.7 minutes**).⁵⁶

Improving endovascular foreign body retrieval with the ÒÑÒ device may reduce costs by avoiding planned and emergency open surgical procedures.

Percutaneous procedures on average incur lower operating time per minute costs than surgical procedures.^{57; 58}

- Mean surgical operating room time cost is approximately **\$24.83 per minute** in the Americas.⁵⁷
- The cost to operate an angiography suite is estimated at **\$539 per hour, or \$8.98 per minute.**⁵⁸



- Surgical ASD closure procedures are estimated to incur higher total costs than transcatheter ASD device closures (2012 USD **\$60,992 vs. \$55,841**). Transcatheter ASD procedures are associated with lower costs for medication, follow-up radiologic and laboratory testing, and a shorter length of stay.⁵⁹
- Index hospital costs are estimated to be significantly lower for transcatheter aortic valve replacements compared to surgical aortic valve replacements, across all patient mortality risk groups (2016–2018; low: **\$61,845 vs. \$68,986**; intermediate: **\$64,658 vs. \$76,965**; high: **\$65,594 vs. \$91,005**).⁶⁰
 - Cost savings were primarily driven by significantly shorter length of stay, which was on average 5 to 8 days shorter across mortality risk groups, for TAVR patients compared to SAVR patients.⁶⁰
- Emergency cardiac procedures can cost approximately **17% to 30% more** than elective procedures.⁶¹

Clinical Evidence

- 1 Coppola J et al. Pediatric Transcatheter Atrial Septal Stent Removal Following ECMO Support Including Use of Novel Retrieval Device. *J Am Coll Cardiol Case Rep.* 2025 Sep, 30 (26). <https://doi.org/10.1016/j.jaccas.2025.104928>
- 2 Elison D, Richard S, Chung C. Percutaneous Removal of a Bullet From the Right Ventricle. *JACC: Cardiovascular Interventions*, 2025, ISSN 1936-8798, <https://doi.org/10.1016/j.jcin.2025.09.002>
- 3 Rago, ARP et al. Successful Use of the ÖNÖ Retrieval Device for Stent Removal in a High-Risk Pediatric Patient. *Journal of the Society for Cardiovascular Angiography & Interventions*, Volume 4, Issue 9, 103812. <https://doi.org/10.1016/j.jscvai.2025.103812>
- 4 Guddeti R et al. Percutaneous Removal of Left Atrial Appendage Occlusion Device With a Dedicated Retrieval System. *J Am Coll Cardiol Case Rep.* 2025 Aug, 30 (24). <https://doi.org/10.1016/j.jaccas.2025.104751>
- 5 Siordia J et al. Retrieval of an Embolized Left Atrial Appendage Occluder From the Aortic Arch Using a Transcatheter Retrieval Device. *J Am Coll Cardiol Case Rep.* 2025 Aug, 30 (24). <https://doi.org/10.1016/j.jaccas.2025.104557>
- 6 Kefer J et al. Management of device embolisation during left atrial appendage closure. *EuroIntervention*, 2025; 21:e838-e846. DOI: 10.4244/EIJ-D-24-00812. <https://eurointervention.pronline.com/article/management-of-device-embolisation-during-left-atrial-appendage-closure>
- 7 Danek B et al. Novel Percutaneous Septectomy for Left Atrial Myxoma. *J Am Coll Cardiol Intv.* 2025 Jul, 18 (13) 1699-1701. <https://doi.org/10.1016/j.jcin.2025.04.030>
- 8 O'Sullivan S et al. Management of Single Leaflet Device Attachment in the Tricuspid Position. *J Am Coll Cardiol Case Rep.* 2025 Jun, 30 (13). <https://doi.org/10.1016/j.jaccas.2025.103864>
- 9 Giustino G et al. First-in-human percutaneous transseptal Retrieval of Embolized Transcatheter valve In the IEft VEtricle (RETRIEVE). *JACC Cardiovasc Interv.* May 06, 2025. Epublished DOI: 10.1016/j.jcin.2025.03.011. <https://doi.org/10.1016/j.jcin.2025.03.011>
- 10 Miller CE, Lim DS, and Shorofsky MJ. Use of an ÖNÖ Retrieval Device for Removal of a Covered CP Stent. *Catheterization and Cardiovascular Interventions*, 2025; 1-2. <https://doi.org/10.1002/ccd.31360>
- 11 Bhatt N et al. Novel Use of a Catheter-Based Device Retrieval System for Percutaneous Treatment of Recurrent Atrial Myxoma. *JACC: Case Reports.* 2025 Jan, 30 (1). <https://www.jacc.org/doi/10.1016/j.jaccas.2024.102737>
- 12 Nienaber S et al. Percutaneous Retrieval of Embolised Left Atrial Appendage Occluder With the Novel ÖNÖ Retrieval Basket. *Canadian Journal of Cardiology.* 2024 Dec; 40(12): 2414 – 2416. [https://onlinecjc.ca/article/S0828-282X\(24\)00928-0/fulltext](https://onlinecjc.ca/article/S0828-282X(24)00928-0/fulltext)
- 13 Kasbaoui S. et al. Percutaneous Excision of Left Atrial Myxoma Using Transcatheter Electrosurgery. *JACC: Case Reports.* Published online December 4, 2024. DOI: 10.1016/j.jaccas.2024.102718. <https://www.jacc.org/doi/10.1016/j.jaccas.2024.102718>
- 14 Chugh Y et al. Percutaneous Retrieval of an Embolized Amulet Device From the Aorta Using a Novel Retrieval System. *JACC Clinical Electrophysiology.* Published online October 2, 2024. DOI: 10.1016/j.jacep.2024.08.007. <https://www.jacc.org/doi/10.1016/j.jacep.2024.08.007>
- 15 Hennawi HA et al. The MICHIGAN Procedure: Mass extraCtIon from the Heart Facilitated by ONOCOR GlobAl embolic protectionN Device. *JACC: Case Reports*, Volume 29, Issue 15, 2024. <https://doi.org/10.1016/j.jaccas.2024.102448>
- 16 Hermann D et al. First-in-Human Percutaneous Removal of Left Atrial Appendage Occlusion Devices With a Novel Retrieval System. *JACC Clinical Electrophysiology.* Published online March 1, 2024. doi: 10.1016/j.jacep.2024.102328. <https://www.jacc.org/doi/10.1016/j.jacep.2024.102328>
- 17 Elison D et al. First-in-Human Percutaneous Excision of a Failed MitraClip Followed by Transcatheter Mitral Valve Replacement. *JACC Cardiovasc Interv.* Published online December 20, 2023. doi:10.1016/j.jcin.2023.11.023. <https://www.jacc.org/doi/10.1016/j.jcin.2023.11.023>
- 18 Taylor AC et al. Retrieval of Large Balloon Fragments During Transcatheter Pulmonary Valve Implantation Using a Novel Retrieval System. *JACC Case Rep.* 2023;26:102058. Published 2023 Oct 18. doi:10.1016/j.jaccas.2023.102058. <https://www.jacc.org/doi/abs/10.1016/j.jaccas.2023.102058>
- 19 Herron C, Batlivala SP, Shahanavaz S. Right atrial thrombus removal with use of the ÖNÖ retrieval device. *Catheter Cardiovasc Interv.* 2023;102(6):1105-1108. doi:10.1002/ccd.30877. <https://onlinelibrary.wiley.com/doi/10.1002/ccd.30877>
- 20 Karunanandaa A et al. Novel use of the ONO Retrieval System for intracardiac thrombectomy in a high-risk pediatric patient. *JSCAI.* 2023 May. <https://doi.org/10.1016/j.jscvai.2023.100976>
- 21 McNamara GPJ et al. Successful removal of a leadless pacemaker from the pulmonary artery via a novel basket retrieval system. *HeartRhythm Case Reports.* 2023 Apr; 9(4): 215-218. <https://doi.org/10.1016/j.hrcr.2022.12.015>
- 22 Steinberg ZL et al. The simplified extraction of atrial tumor with targeted loop electricity (SEATTLE) procedure. *JACC: Case Reports.* 2023 Mar. <https://www.jacc.org/doi/10.1016/j.jaccas.2023.101758>

For complete list of clinical evidence, visit our website.



References

- 1 B Braun Interventional Systems Inc (2023). ONO Retrieval Device. B Braun Interventional Systems Inc: Available online at: <https://www.bisusa.com/en/products/b/ono-retrieval-device.html>. Accessed: 2024/02/14.
- 2 ONOCOR (Web Page) ÕNÛCOR 12F ÕNÛ RETRIEVAL DEVICE: Instructions for Use. Updated Available online at: <https://www.onocorvascular.com/instructions-for-use#:~:text=The%20Onocor%20%C5%8C%20Retrieval%20device,using%20minimally%20invasive%20surgical%20procedures>. Accessed: 2024/02/14.
- 3 Onocor (2022). Data on File: Onocor Value and Pricing Assessment Final report. Unpublished
- 4 Correia R, Garcia A, Camacho N, Catarino J, Bento R et al. (2022) Intravascular foreign body retrieval. *Port J Card Thorac Vasc Surg* 29 (2): 45-50.
- 5 Lyu T, Cao S, Wang J, Song L, Tong X et al. (2022) Endovascular Removal of Foreign Bodies: Single Center Experience. *Indian Journal of Surgery* 85 (Suppl 1): 254-258.
- 6 Windecker S, Gilard M, Achenbach S, Cribier A, Delgado V et al. (2024) Device innovation in cardiovascular medicine: a report from the European Society of Cardiology Cardiovascular Round Table. *Eur Heart J* 45 (13): 1104-1115.
- 7 Alkhouli M, Sievert H, Rihal CS (2019) Device Embolization in Structural Heart Interventions: Incidence, Outcomes, and Retrieval Techniques. *JACC Cardiovasc Interv* 12 (2): 113-126.
- 8 Schechter MA, O'Brien PJ, Cox MW (2013) Retrieval of iatrogenic intravascular foreign bodies. *J Vasc Surg* 57 (1): 276-281.
- 9 Leite TFO, Pazinato LV, Bortolini E, Pereira OI, Nomura CH et al. (2022) Endovascular Removal of Intravascular Foreign Bodies: A Single-Center Experience and Literature Review. *Ann Vasc Surg* 82 362-376.
- 10 Surov A, Wienke A, Carter JM, Stoevesandt D, Behrmann C et al. (2009) Intravascular embolization of venous catheter--causes, clinical signs, and management: a systematic review. *JPEN J Parenter Enteral Nutr* 33 (6): 677-685.
- 11 McNamara GPJ, Haber ZM, Lee EW, Levi DS, Bender A (2023) Successful removal of a leadless pacemaker from the pulmonary artery via a novel basket retrieval system. *HeartRhythm Case Rep* 9 (4): 215-218.
- 12 Mack M, Carroll JD, Thourani V, Vemulapalli S, Squiers J et al. (2022) Transcatheter Mitral Valve Therapy in the United States: A Report from the STS/ACC TVT Registry. *Ann Thorac Surg* 113 (1): 337-365.
- 13 Chhatriwalla AK, Vemulapalli S, Holmes DR, Jr, Dai D, Li Z et al. (2019) Institutional Experience With Transcatheter Mitral Valve Repair and Clinical Outcomes: Insights From the TVT Registry. *JACC Cardiovasc Interv* 12 (14): 1342-1352.
- 14 Vishwanath V, Akseer S, Frankfurter C, Reza S, Farrell A et al. (2022) Comparative effectiveness of devices for transcatheter closure of atrial septal defects: Systematic review and network meta-analysis. *Arch Cardiovasc Dis* 115 (12): 664-674.
- 15 Garre S, Gadhinglajkar S, Sreedhar R, Krishnamoorthy KM, Pillai VV (2023) Atrial septal defect occluder device embolization: Experience of a tertiary care cardiac center. *Ann Card Anaesth* 26 (2): 149-154.
- 16 Friedman DJ, Freeman JV, Zimmerman S, Tan Z, Pereira L et al. (2023) Watchman device migration and embolization: A report from the NCDR LAAO Registry. *J Cardiovasc Electrophysiol* 34 (5): 1192-1195.
- 17 Price MJ, Friedman DJ, Du C, Wang Y, Lin Z et al. (2022) Comparative Safety of Transcatheter LAAO With the First-Generation Watchman and Next-Generation Watchman FLX Devices. *JACC Cardiovasc Interv* 15 (21): 2115-2123.
- 18 Badesha AS, Siddiqui MM, Bains BRS, Bains PRS, Khan T (2022) A Systematic Review on the Incidence of Stent Migration in the Treatment of Acute and Chronic Iliofemoral Disease Using Dedicated Venous Stents. *Ann Vasc Surg* 83 328-348.
- 19 Seager MJ, Busuttill A, Dharmarajah B, Davies AH (2016) Editor's Choice-- A Systematic Review of Endovenous Stenting in Chronic Venous Disease Secondary to Iliac Vein Obstruction. *Eur J Vasc Endovasc Surg* 51 (1): 100-120.
- 20 Palmisano P, Iacopino S, De Vivo S, D'Agostino C, Tomasi L et al. (2022) Leadless transcatheter pacemaker: Indications, implantation technique and peri-procedural patient management in the Italian clinical practice. *Int J Cardiol* 365 49-56.
- 21 Schiavone M, Filtz A, Gasperetti A, Breitenstein A, Palmisano P et al. (2022) Leadless Pacemaker Implantation in the Emergency Bradycardia Setting: Results from a Multicenter European Registry. *Medicina (Kaunas)* 59 (1): 22
- 22 Darlington D, Brown P, Carvalho V, Bourne H, Mayer J et al. (2022) Efficacy and safety of leadless pacemaker: A systematic review, pooled analysis and meta-analysis. *Indian Pacing Electrophysiol J* 22 (2): 77-86.
- 23 Roberts PR, Clementy N, Al Samadi F, Garweg C, Martinez-Sande JL et al. (2017) A leadless pacemaker in the real-world setting: The Micra Transcatheter Pacing System Post-Approval Registry. *Heart Rhythm* 14 (9): 1375-1379.
- 24 Ali S, Lee SR, Strosberg D, Aboian E, Guzman R et al. (2023) Balloon rupture with eversion during innominate vein angioplasty requiring surgical retrieval. *J Vasc Surg Cases Innov Tech* 9 (3): 101242.
- 25 Seropian IM, Romeo FJ, Falconi M, Agatiello CR, Berrocal DH (2020) Balloon rupture during aortic valvuloplasty with compliant balloon: predictors and outcomes. *Cardiovasc Interv Ther* 35 (3): 291-299.
- 26 Herron C, Batlivala SP, Shahanavaz S (2023) Right atrial thrombus removal with use of the ONO retrieval device. *Catheter Cardiovasc Interv* 102 (6): 1105-1108.
- 27 Karunanandaa A, Patel A, Wong P, Vo C, Van Tienderen M et al. (2023) Novel Use of the ÕNÛ Retrieval System for Intracardiac Thrombectomy in a High-Risk Pediatric Patient. *Journal of the Society for Cardiovascular Angiography & Interventions* 2 (4): 100976.
- 28 Steinberg ZL, Elison D, Vincent LL, Oxorn D, McCabe JM (2023) The Simplified Extraction of Atrial Tumor With Targeted Loop Electricity (SEATTLE) Procedure. *JACC Case Rep* 10 101758.
- 29 Shanmugasundaram M, Chatterjee A, Lee K (2022) Percutaneous Large Thrombus and Vegetation Evacuation in the Catheterization Laboratory. *Interv Cardiol Clin* 11 (3): 349-358.
- 30 Keeton J, Gonzalez PE, Cox J, Morlend RM, Amin AA et al. (2020) Saved by the VAC: Minimally Invasive Removal of a Surely Fatal Right Heart Thrombus in a Patient with Advanced Heart Failure. *Case Rep Cardiol* 2020 7579262.
- 31 Athappan G, Sengodan P, Chacko P, Gandhi S (2015) Comparative efficacy of different modalities for treatment of right heart thrombi in transit: a pooled analysis. *Vasc Med* 20 (2): 131-138.
- 32 Jolugbo P, Ariens RAS (2021) Thrombus Composition and Efficacy of Thrombolysis and Thrombectomy in Acute Ischemic Stroke. *Stroke* 52 (3): 1131-1142.
- 33 Czapllicki C, Albadawi H, Partovi S, Gandhi RT, Quencer K et al. (2017) Can thrombus age guide thrombolytic therapy? *Cardiovasc Diagn Ther* 7 (Suppl 3): S186-S196.
- 34 Tran MH, Wilcox T, Tran PN (2020) Catheter-related right atrial thrombosis. *J Vasc Access* 21 (3): 300-307.
- 35 Tyebally S, Chen D, Bhattacharyya S, Mughrabi A, Hussain Z et al. (2020) Cardiac Tumors: JACC CardioOncology State-of-the-Art Review. *JACC CardioOncol* 2 (2): 293-311.
- 36 Rahoma M, Arisha MJ, Elmously A, El-Sayed Ahmed MM, Spadaccio C et al. (2020) Cardiac tumors prevalence and mortality: A systematic review and meta-analysis. *Int J Surg* 76 178-189.
- 37 Moscarelli M, Rahoma M, Nasso G, di Bari N, Speciale G et al. (2021) Minimally invasive approaches to primary cardiac tumors: A systematic review and meta-analysis. *J Card Surg* 36 (2): 483-492.
- 38 Beurskens NEG, Breeman KTN, Dasseljaar KJ, Meijer AC, Quast ABE et al. (2019) Leadless cardiac pacing systems: current status and future prospects. *Expert Rev Med Devices* 16 (11): 923-930.
- 39 Beurskens NE, Tjong FV, Knops RE (2017) End-of-life Management of Leadless Cardiac Pacemaker Therapy. *Arrhythm Electrophysiol Rev* 6 (3): 129-133.
- 40 Li J, Hou WB, Cao MK, Zhou WX, Wang Y et al. (2019) Safety and efficacy of leadless pacemaker retrieval. *J Cardiovasc Electrophysiol* 30 (9): 1671-1678.
- 41 B Braun Interventional Systems Inc (2023). Data on File: ÕNÛ™ Retrieval Device Unpublished Slide Presentation.
- 42 Koussayer S, Mousa A, Elkalla MA, Alaujan R, Khalil B (2021) Conservative management of asymptomatic lately embolized atrial septal occluder device to the supraceliac abdominal aorta: case report and the literature review. *Radiol Case Rep* 16 (7): 1912-1918.
- 43 Ooi YK, Kelleman M, Ehrlich A, Glanville M, Porter A et al. (2016) Transcatheter Versus Surgical Closure of Atrial Septal Defects in Children: A Value Comparison. *JACC Cardiovasc Interv* 9 (1): 79-86.
- 44 Felbel D, Paukovitsch M, Forg R, Stephan T, Mayer B et al. (2022) Comparison of transcatheter edge-to-edge and surgical repair in patients with functional mitral regurgitation using a meta-analytic approach. *Front Cardiovasc Med* 9 1063070.
- 45 Sarris GE, Kirvasilis G, Zavaropoulos P, Belli E, Berggren H et al. (2010) Surgery for complications of trans-catheter closure of atrial septal defects: a multi-institutional study from the European Congenital Heart Surgeons Association. *Eur J Cardiothorac Surg* 37 (6): 1285-1290.
- 46 Conzelmann L, Grotherr P, Dapeng L, Wurth A, Widder J et al. (2023) Valve-Related Complications in TAVI Leading to Emergent Cardiac Surgery. *Thorac Cardiovasc Surg* 71 (2): 107-117.
- 47 Khan MR, Kayani WT, Manan M, Munir A, Hamzeh I et al. (2020) Comparison of surgical versus transcatheter aortic valve replacement for patients with aortic stenosis at low-intermediate risk. *Cardiovasc Diagn Ther* 10 (2): 135-144.
- 48 B Braun Interventional Systems Inc (2023). Data on File: ONO Retrieval Device Clinical Cases. Unpublished
- 49 Hermann D, Khan Z, Khan Z, Abreu J, Loyalka P et al. (2024) First in Man Percutaneous Removal of Left Atrial Appendage Occlusion Devices with a Novel Retrieval System. *JACC: Clinical Electrophysiology*
- 50 Taylor AC, Ghandour MAH, Khan A, Gowda ST, Nunez-Gallegos F et al. (2023) Retrieval of Large Balloon Fragments During Transcatheter Pulmonary Valve Implantation Using a Novel Retrieval System. *JACC Case Rep* 26 102058.
- 51 Elison D, Aldea G, Jelacic S, Chung CJ, Mackensen GB et al. (2023) First-in-Human Percutaneous Excision of a Failed MitraClip Followed by Transcatheter Mitral Valve Replacement. *JACC Cardiovasc Interv*
- 52 Rader F, Siegel RJ, Flint N, Minhas H, Shmueli YG et al. (2021) Mitral Valve Surgery After Failed MitraClip: A Single-Center Experience. *J Invasive Cardiol* 33 (4): E236-E243.
- 53 Mkalaluh S, Szczechowicz M, Torabi S, Schmack B, Sabashnikov A et al. (2017) Surgical Treatment of Cardiac Tumors: Insights from an 18-Year Single-Center Analysis. *Med Sci Monit* 23 6201-6209.
- 54 Pazinato LV, Leite TFO, Bortolini E, Pereira OI, Nomura CH et al. (2022) Percutaneous retrieval of intravascular foreign body in children: a case series and review. *Acta Radiol* 63 (5): 684-691.
- 55 Sundermann SH, Sromicki J, Rodriguez Cetina Bieffer H, Seifert B, Holubec T et al. (2014) Mitral valve surgery: right lateral minithoracotomy or sternotomy? A systematic review and meta-analysis. *J Thorac Cardiovasc Surg* 148 (5): 1989-1995 e1984.
- 56 Sondergaard L, Popma JJ, Reardon MJ, Van Mieghem NM, Deeb GM et al. (2019) Comparison of a Complete Percutaneous Versus Surgical Approach to Aortic Valve Replacement and Revascularization in Patients at Intermediate Surgical Risk: Results From the Randomized SURTAVI Trial. *Circulation* 140 (16): 1296-1305.
- 57 Chen BP, Clymer JW, Turner AP, Ferko N (2019) Global hospital and operative costs associated with various ventral cavity procedures: a comprehensive literature review and analysis across regions. *J Med Econ* 22 (11): 1210-1220.
- 58 Beheshti MV, Meek J (2014) Calculation of operating expenses for conventional transarterial chemoembolization in an academic medical center: a step toward defining the value of transarterial chemoembolization. *J Vasc Interv Radiol* 25 (4): 567-574.
- 59 O'Byrne ML, Gillespie MJ, Shinohara RT, Dori Y, Rome JJ et al. (2015) Cost comparison of transcatheter and operative closures of ostium secundum atrial septal defects. *Am Heart J* 169 (5): 727-735 e722.
- 60 Baron SJ, Ryan MP, Moore KA, Clancy SJ, Gunnarsson CL (2022) Contemporary Costs Associated With Transcatheter Versus Surgical Aortic Valve Replacement in Medicare Beneficiaries. *Circ Cardiovasc Interv* 15 (3): e011295.
- 61 Haider AH, Obirieze A, Velopulos CG, Richard P, Latif A et al. (2015) Incremental Cost of Emergency Versus Elective Surgery. *Ann Surg* 262 (2): 260-266.

The ÕNÕ Retrieval Device is indicated for use in the cardiovascular system to retrieve foreign objects using minimally invasive surgical procedures. Refer to the Instructions for Use for complete indications, relevant warnings, precautions, complications and contraindications. ÕNÕ and ÕNÕCOR are registered trademarks of ÕNÕCOR LLC.

For more information or to place an order, contact your B. Braun Interventional Systems Inc. Representative or call 877-836-2228.
Rx only
Distributed by: B. Braun Interventional Systems Inc. | a member of the B. Braun Group of Companies
Bethlehem, PA 18018 | USA | Tel 877-836-2228 | Fax 610-849-1334 | www.bisusa.com

