# B. Braun Inflation Device

## GB Instructions for Use

**INDICATION**  
This Inflation Device is recommended for use during the performance of balloon angioplasty procedures to inflate and deflate the balloon and to measure the pressure within the balloon during the procedure.

**DESCRIPTION**  
The disposable B. Braun Inflation Device consists of a 25 cc syringe barrel with a pressure gauge, a proximal pressure assembly with Winged Lock Mechanism used to generate and control balloon inflation pressures, a flexible high pressure extension tube with proximal pressure stopcock, the gauge is calibrated from 0 to 30 atm (0 to 441 psi) of positive pressure. The accuracy of the gauge is within 1 atm over the range. Accuracy of syringe graduations: ±5%.

**WARNINGS**

- Isowave 370 or Renovist® II in combination with saline may cause renal complications. Patients should be used with Renovist® patients. Refer to the directions accompanying the balloon dilation catheter for instructions for use, maximum pressure, precautions, and warnings for that device.
- Re-use of single-use devices carries a potential risk of patient or user. It may lead to contamination and/or impairment of functional capacity. Contamination and/or limited functionality of the device may lead to injury. Inflex or death of the patient.

**CAUTION:** Read the instructions for use of all equipment to be used in the procedure prior to use.

**PRECAUTIONS**

- Follow all instructions recommended by the balloon dilation catheter manufacturer.
- Use the sterile contrast media that the balloon dilation catheter manufacturer recommends.
- The balloon dilation catheter is a single use device. Its recommended maximum balloon inflation pressure should not be exceeded.
- Do NOT exceed a pressure of 441 psi (30 atm) when using this device.

**WARNING:** Refer to the manufacturer's directions accompanying the balloon dilation catheter for specific information on use, maximum inflation pressure, precautions, and warnings for that device. Inflation pressures should be closely monitored when inflating balloons. The inflation syringe is a high volume, low compliance system capable of generating high pressures with relative ease.

**INSTRUCTIONS FOR USE:** Use Aseptic Technique.

- Inspect prior to use to ensure no damage has occurred during shipping and handling. Do not use if packaging or product is damaged or contains any of protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness or death of the patient.

**NOTE:** Product is packaged in a sterile prep pouch to facilitate sterile transfer of the device. Inspect for damage prior to use. Follow steps below to check inflation device.

- Dispose of according to local guidelines and/or clinical protocol.
- If during the use of this product as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.
- The winged lock mechanism counter-clockwise so wings are parallel to the gauge.

2. Aspirate air by pulling plunger back to 25 cc graduation, 3. Attach stopcock to the proximal end of inflation device and turn stopcock handle with the "WINGED" position facing extension tubing port.

4. Lock syringe by turning the winged lock mechanism clockwise 90° so wings are perpendicular to the gauge.

5. Turn the plunger handle clockwise and advance the syringe plunger forward to build pressure.

6. Measure inflation pressure to approximately 441 psi (30 atm), let gauge stabilize 10 seconds and observe gauge for any loss of pressure.

7. Do not use if damage or loss of pressure is observed.

**PREPARATION**

1. Prepare solution of contrast media and normal saline as recommended by catheter manufacturer.
2. Insert the end of the extension tube into the contrast solution.
3. Aspirate up to 20 cc of solution into the Inflation Syringe.
4. Purge all air from the inflation syringe barrel by pointing it upwards so that the air accumulates in the distal end of the barrel. Slowly advance the syringe plunger until all air has been expelled from the syringe barrel and extension tube.

**CAUTION: INSPECT THE SYRINGE TUBING AND STOPCOCK (IF USED) TO INSURE THAT THERE IS NO AIR IN THE SYSTEM, RECOMMENDED OPTIMAL VOLUME OF CONTRAST SOLUTION IS 10 TO 10 cc.**

**ATTACHING THE INFLATION DEVICE TO THE BALLOON**

1. Prepare and test the balloon catheter according to the manufacturer's instructions for use.
2. Create a fluid-tight connection between the balloon and the Inflation Syringe extension tube and connect the luer connectors securely.
3. Unlock by turning the Winged Lock Mechanism 90° counter-clockwise and lock back on the plunger handle to apply a vacuum to the balloon.

**INFLATION INSTRUCTIONS**

1. To engage, turn the winged lock mechanism clockwise 90°.
2. Ensure the latch is fully engaged.
3. Slowly rotate the plunger clockwise to increase the pressure of contrast media to the desired pressure.

**WARNING: A SLIGHT ROTATION OF THE HANDLE WILL CHANGE THE PRESSURE SIGNIFICANTLY. MONITOR PRESSURES CONTINUALLY TO AVOID OVER PRESSURIZATION.**

**DEFLATION INSTRUCTIONS**

Release the winged lock mechanism counter-clockwise 90° to a horizontal position.

Rotate the plunger fully withdrawn and all fluid is removed from the balloon.

**NOTE:** Negative pressure can be maintained by engaging the winged lock mechanism when the plunger is fully withdrawn.

**CAUTION: OBSERVE THE ANALOG GAUGE TO CONFIRM THAT THE SYSTEM MAINTAINS A VACUUM. LOSS OF NEGATIVE PRESSURE INDICATES THE PRESENCE OF A LEAK IN THE SYSTEM.**

Isowave is a registered trademark of Bracco Diagnostics, Inc. Renovist is a trademark of Bracco Diagnostics, Inc.

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B. Braun Medical Inc.
B. Braun Melsungen AG