SAFECAN™ Safety Introducer Needle

CE 0123

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GB **Instructions for Use**
Indications
The Safety Introducer Needle incorporates an active safety clip mechanism to minimize inadvertent needlesticks when the needle is used to introduce a guidewire.
Contraindications
The Safety Introducer Needle should not be used in patients with known hypersensitivity to any of the materials employed.
Potential Complications
Potential complications with the Safety Introducer Needle include, but are not limited to, the following:
• Inadvertent vessel puncture • Pneumothorax
• Hemothorax • Air embolism
• Hematoma formation • Brachial plexus injury
• Thrombosis
Risks
This Safety Introducer Needle is designed to reduce the risk of accidental needlesticks; however, care must be taken to avoid needlesticks. Universal precautions must be adhered to in accordance with Centers for Disease Control and Prevention/Occupational Safety and Health Administration (CDC/OSHA) standards for bloodborne pathogens, when utilizing any needles that penetrates the skin to avoid the risk of exposure to contaminated blood.

WARNINGS
• After withdrawal of the Safety Introducer Needle from the patient and activation of the safety clip, do not attempt to re-insert the needle into the patient. Also, do not attempt to recap the needle with the protective guard. Dispose of the needle immediately into an appropriate sharps container.
• Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
Instructions for Use: Use Aseptic Technique.
• **NOTE:** Both hands should remain behind the needle bevel at all times.
• Inspect prior to use to insure no damage has occurred during shipping and handling. Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness or death of the patient.
• If during the use of this product or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.
1. After disinfection of the puncture site and removal of the protective guard from the needle, puncture desired vessel.
2. Insert guidewire into the proximal end of the needle hub.
3. Advance guidewire through needle into vessel.
4. Advancement of guidewire may require gentle rotating motion.
CAUTION: DO NOT WITHDRAW GUIDEWIRE AGAINST NEEDLE BEVEL TO AVOID POSSIBLE SEVERING OR DAMAGING OF GUIDEWIRE. IF GUIDEWIRE WILL NOT ADVANCE, SIMULTANEOUSLY REMOVE NEEDLE AND GUIDEWIRE FROM PUNCTURE SITE AND RESTART.

5. Advance guidewire in routine fashion to desired depth.
6. Hold guidewire in place, withdraw Safety Introducer Needle from the vessel and remove the needle from guidewire.
CAUTION: DO NOT ACTIVATE SAFETY CLIP BEFORE REMOVING THE NEEDLE FROM GUIDEWIRE.
7. Grasp clear clip housing, which houses the metal safety clip, with your free hand. Advance clear clip housing forward along the Safety Introducer Needle toward needle bevel.
8. In one continuous motion, advance clip housing over needle bevel activating safety mechanism. The clip will remain on the needle bevel. The clip housing will separate from the clip and needle.
9. Dispose of needle and clear clip housing immediately into sharps container.

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