Celsite®

Implantable Access Port System
(Venous)

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Table of Contents

Glossary.....	1
Purpose of the device	2
Description of the device	3
How will the port be implanted?.....	4
Expectations of the device	5
When the device should not be accessed.....	6
Benefits and risks of the device.....	6
Patient Identification Card	7

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Glossary

Your physician may use these terms when describing your treatment or care of the Celsite® Implantable Access Port System.

Catheter: A thin plastic tube that delivers medication from the port.

Access Port: The triangular shaped device with a septum where the needle is inserted.

Heparin: A medication that helps keep the implanted port and catheter free of blood clots.

Septum: Made from a self sealing material which allows the needle to access the port.


Port site: Area of the body where the port is implanted.

Purpose of the device

The Celsite® Implantable Access Port System (Celsite® Access Port) is a small titanium and plastic implantable port designed to allow safe, repeated access for administration of chemotherapy, parenteral nutrition, blood sampling or transfusions, and administration of antibiotics or antiviral drugs.


The Celsite® Access Port is completely covered by your skin and offers you total freedom of movement.

This booklet is provided to help you understand how the Celsite® Access Port will be used in your treatment. This information is provided to you as a reference only. We urge you to discuss any questions you have regarding the device or treatment with your doctor or nurse.

Description of the device

The Celsite® Access Port has two components: a chamber called a **port**, and a thin plastic tube called a **catheter**. The port is a small triangular shaped device which has a raised center referred to as the **septum**. The septum is made from a self-sealing material. The septum is where the needle is inserted for delivery of medication and fluids. The medication is carried from the port through the catheter. All ports and catheters are latex and DEHP free.

1. Silicone septum

2. Plastic or Epoxy resin covering

3. Titanium chamber


4. Silicone or polyurethane catheter

How will the port be implanted?

A brief outpatient procedure is performed usually under local anesthesia. This means you will most likely be able to return home on the same day. The doctor will place the port just beneath the skin in a position most suitable for your treatment. For many patients, the port is placed on the upper part of the chest. The end of the catheter will be placed in a large blood vessel for venous access. With the system in place, you will see only a small bump under your skin.

After your incision(s) heals, your doctor will let you know when you may resume your normal activities. Your skin's natural barrier reduces the chance of infection or damage to the port.

Call your doctor before visiting the dentist and ask if you should get a prescription for antibiotics beginning the day before and ending the day after dental work.

Possible sites where a port may be placed

Expectations of the device

Your doctor or nurse can insert a needle through the skin and into the septum of the port. You may feel a slight prick as the needle goes through your skin. Topical anesthesia may be used at your physician discretion. This needle can then be used to deliver medications or other fluids, as well as withdrawal of blood samples for laboratory testing.

Injections of medications can be done over a very short period of time (called a push or bolus injection) or over a prolonged period of time (called continuous infusion). A continuous infusion may last anywhere from one hour up to several days or longer. A dressing will usually be placed over the needle to stabilize it in the port and a small tube will connect the needle to the medication source.

After each treatment, and occasionally between treatments, your doctor or nurse will rinse the port with a small amount of **heparin** solution to help keep the catheter clear of blood clots. If your treatments are spaced more than a month apart, you may need to have the port rinsed between treatments by your nurse. Your doctor will give you a schedule for these rinses.

When the device should not be accessed

It is important to examine the skin around the Celsite® Access Port daily to detect any changes in the appearance of the port implantation site. Remember: some redness and tenderness is normal in the area of the Celsite® Access Port for several days following implantation.

Contact your doctor or nurse **immediately** if you notice any of the following symptoms of infection or other complications related to the port implantation.

- **Fever, chills, or other flu-like symptoms**
- **Swelling, pain, pus, redness or burning sensation around the port incision**
- **Shortness of breath or dizziness**
- **Chest pain**
- **Pain or swelling of neck, face, or arm on side of body where the implantable port is inserted**

Resist the urge to manipulate the port or the skin around the port. This could break the small stitches that hold the port in place or cause an infection.

Benefits and risks of the device

Your Celsite® Access Port provides a means for your clinician to easily deliver the medications and fluids that you need over an extended period of time.

The access port, because it is implanted under your skin, is more comfortable and may carry a reduced risk of infection compared to external catheters. Patients with a Celsite® Access Port may safely resume activities such as swimming, bathing, jogging etc. after the incision heals.

Potential risks associated with implantable Access Port Systems include, but are not limited to infection, catheter blockage or breakage, and swelling or pain in the area around the system.

Patient Identification Card

After the implantation of your access port you should have received a patient identification card. It is very important for people who have implanted devices to always carry the port patient identification card with pertinent medical information. This is particularly important in the event of an emergency, so medical personnel will be aware you have a port.

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en - Celsite® is a Registered Trademark.
B. Braun Interventional Systems Inc.