

BIB® CP Stent® Placement Catheter

Reference Number	Outer Balloon Diameter (mm)	Outer Balloon Length (cm)	Outer Rated Burst Pressure (atm)	Inner Balloon Diameter (mm)	Inner Balloon Length (cm)	Inner Rated Burst Pressure (atm)	Intro Size (Fr)	Shaft Size (Fr)
614402	12.0	2.5	7.0	6.0	1.5	5.0	8.0	8.0
614403	12.0	3.0	7.0	6.0	2.0	5.0	8.0	8.0
614404	12.0	3.5	7.0	6.0	2.5	5.0	8.0	8.0
614405	12.0	4.0	7.0	6.0	3.0	5.0	8.0	8.0
614406	12.0	4.5	7.0	6.0	3.5	5.0	8.0	8.0
614407	12.0	5.0	7.0	6.0	4.0	5.0	8.0	8.0
614408	12.0	5.5	7.0	6.0	4.5	5.0	8.0	8.0
614409	14.0	2.5	6.0	7.0	1.5	5.0	8.0	8.0
614410	14.0	3.0	6.0	7.0	2.0	5.0	8.0	8.0
614411	14.0	3.5	6.0	7.0	2.5	5.0	8.0	8.0
614412	14.0	4.0	6.0	7.0	3.0	5.0	8.0	8.0
614413	14.0	4.5	6.0	7.0	3.5	5.0	8.0	8.0
614414	14.0	5.0	6.0	7.0	4.0	5.0	8.0	8.0
614415	14.0	5.5	6.0	7.0	4.5	5.0	8.0	8.0
614416	15.0	2.5	5.0	7.0	1.5	5.0	9.0	9.0
614417	15.0	3.0	5.0	7.0	2.0	5.0	9.0	9.0
614418	15.0	3.5	5.0	7.0	2.5	5.0	9.0	9.0
614419	15.0	4.0	5.0	7.0	3.0	5.0	9.0	9.0
614458	15.0	4.5	5.0	7.0	3.5	5.0	9.0	9.0
614421	15.0	5.0	5.0	7.0	4.0	5.0	9.0	9.0
614422	15.0	5.5	5.0	7.0	4.5	5.0	9.0	9.0
614423	16.0	2.5	5.0	8.0	1.5	5.0	9.0	9.0
614424	16.0	3.0	5.0	8.0	2.0	5.0	9.0	9.0
614425	16.0	3.5	5.0	8.0	2.5	5.0	9.0	9.0
614426	16.0	4.0	5.0	8.0	3.0	5.0	9.0	9.0
614427	16.0	4.5	5.0	8.0	3.5	5.0	9.0	9.0
614428	16.0	5.0	5.0	8.0	4.0	5.0	9.0	9.0
614429	16.0	5.5	5.0	8.0	4.5	5.0	9.0	9.0
614627	16.0	6.0	5.0	8.0	5.0	5.0	9.0	9.0
614628	16.0	6.5	5.0	8.0	5.5	5.0	9.0	9.0
614430	18.0	2.5	4.0	9.0	1.5	5.0	10.0	9.0
614431	18.0	3.0	4.0	9.0	2.0	5.0	10.0	9.0
614432	18.0	3.5	4.0	9.0	2.5	5.0	10.0	9.0
614433	18.0	4.0	4.0	9.0	3.0	5.0	10.0	9.0
614434	18.0	4.5	4.0	9.0	3.5	5.0	10.0	9.0
614435	18.0	5.0	4.0	9.0	4.0	5.0	10.0	9.0
614436	18.0	5.5	4.0	9.0	4.5	5.0	10.0	9.0
614629	18.0	6.0	4.0	9.0	5.0	5.0	10.0	9.0
614630	18.0	6.5	4.0	9.0	5.5	5.0	10.0	9.0
614437	20.0	3.0	4.0	10.0	2.0	5.0	10.0	9.0
614438	20.0	3.5	4.0	10.0	2.5	5.0	10.0	9.0
614439	20.0	4.0	4.0	10.0	3.0	5.0	10.0	9.0
614440	20.0	4.5	4.0	10.0	3.5	5.0	10.0	9.0
614441	20.0	5.0	4.0	10.0	4.0	5.0	10.0	9.0
614442	20.0	5.5	4.0	10.0	4.5	5.0	10.0	9.0
614631	20.0	6.0	4.0	10.0	5.0	5.0	10.0	9.0
614632	20.0	6.5	4.0	10.0	5.5	5.0	10.0	9.0
614443	22.0	3.0	3.0	11.0	2.0	4.5	11.0	9.0
614444	22.0	3.5	3.0	11.0	2.5	4.5	11.0	9.0
614445	22.0	4.0	3.0	11.0	3.0	4.5	11.0	9.0
614446	22.0	4.5	3.0	11.0	3.5	4.5	11.0	9.0
614447	22.0	5.0	3.0	11.0	4.0	4.5	11.0	9.0
614448	22.0	5.5	3.0	11.0	4.5	4.5	11.0	9.0
614633	22.0	6.0	3.0	11.0	5.0	4.5	11.0	9.0
614634	22.0	6.5	3.0	11.0	5.5	4.5	11.0	9.0
614449	24.0	3.0	3.0	12.0	2.0	4.5	11.0	9.0
614450	24.0	3.5	3.0	12.0	2.5	4.5	11.0	9.0
614451	24.0	4.0	3.0	12.0	3.0	4.5	11.0	9.0
614452	24.0	4.5	3.0	12.0	3.5	4.5	11.0	9.0
614453	24.0	5.0	3.0	12.0	4.0	4.5	11.0	9.0
614454	24.0	5.5	3.0	12.0	4.5	4.5	11.0	9.0
614469	24.0	6.0	3.0	12.0	5.0	4.5	11.0	9.0
614635	24.0	6.5	3.0	12.0	5.5	4.5	11.0	9.0
614459	26.0	4.0	3.0	13.0	3.0	4.0	16.0	11.0
614460	26.0	5.0	3.0	13.0	4.0	4.0	16.0	11.0
614461	26.0	6.0	3.0	13.0	5.0	4.0	16.0	11.0
614636	26.0	6.5	3.0	13.0	5.5	4.0	16.0	11.0
614462	28.0	4.0	2.0	14.0	3.0	4.0	16.0	11.0
614463	28.0	5.0	2.0	14.0	4.0	4.0	16.0	11.0
614464	28.0	6.0	2.0	14.0	5.0	4.0	16.0	11.0
614637	28.0	6.5	2.0	14.0	5.5	4.0	16.0	11.0
614465	30.0	4.0	2.0	15.0	3.0	4.0	16.0	11.0
614466	30.0	5.0	2.0	15.0	4.0	4.0	16.0	11.0
614467	30.0	6.0	2.0	15.0	5.0	4.0	16.0	11.0
614638	30.0	6.5	2.0	15.0	5.5	4.0	16.0	11.0

All BIB catheters have a usable length of 110cm and 0.035" guidewire compatibility.

BIB (BALLOON IN BALLOON) CATHETERS - ORDERING INFORMATION



BIB® CP Stent® Placement Catheter

Reference Number	Outer Balloon Diameter (mm)	Outer Balloon Length (cm)	Outer Rated Burst Pressure (atm)	Inner Balloon Diameter (mm)	Inner Balloon Length (cm)	Inner Rated Burst Pressure (atm)	Intro Size (Fr)	Shaft Size (Fr)
614621	12.0	6.0	7.0	6.0	5.0	5.0	8.0	8.0
614622	12.0	6.5	7.0	6.0	5.5	5.0	8.0	8.0
614623	14.0	6.0	6.0	7.0	5.0	5.0	8.0	8.0
614624	14.0	6.5	6.0	7.0	5.5	5.0	8.0	8.0
614625	15.0	6.0	5.0	7.0	5.0	5.0	9.0	9.0
614626	15.0	6.5	5.0	7.0	5.5	5.0	9.0	9.0

All BIB catheters have a usable length of 110cm and 0.035" guidewire compatibility.

BIB® PTA Dilatation Catheter

Reference Number	Outer Balloon Diameter (mm)	Outer Balloon Length (cm)	Outer Rated Burst Pressure (atm)	Inner Balloon Diameter (mm)	Inner Balloon Length (cm)	Inner Rated Burst Pressure (atm)	Intro Size (Fr)	Shaft Size (Fr)
614388	8.0	2.5	10.0	4.0	1.5	5.0	8.0	8.0
614389	8.0	3.0	10.0	4.0	2.0	5.0	8.0	8.0
614390	8.0	3.5	10.0	4.0	2.5	5.0	8.0	8.0
614391	8.0	4.0	10.0	4.0	3.0	5.0	8.0	8.0
614392	8.0	4.5	10.0	4.0	3.5	5.0	8.0	8.0
614393	8.0	5.0	10.0	4.0	4.0	5.0	8.0	8.0
614394	8.0	5.5	10.0	4.0	4.5	5.0	8.0	8.0
614395	10.0	2.5	9.0	5.0	1.5	5.0	8.0	8.0
614396	10.0	3.0	9.0	5.0	2.0	5.0	8.0	8.0
614397	10.0	3.5	9.0	5.0	2.5	5.0	8.0	8.0
614398	10.0	4.0	9.0	5.0	3.0	5.0	8.0	8.0
614399	10.0	4.5	9.0	5.0	3.5	5.0	8.0	8.0
614455	10.0	5.0	9.0	5.0	4.0	5.0	8.0	8.0
614456	10.0	5.5	9.0	5.0	4.5	5.0	8.0	8.0

All BIB catheters have a usable length of 110cm and 0.035" guidewire compatibility.

***Lead times are estimated.**

Please contact your local B. Braun Interventional Systems Medical Device Specialist or Customer Support Representative to confirm availability and delivery.

Indications for Use:

The BIB® CP Stent® Placement Catheter is indicated for CP Stent®/Covered CP Stent™ placement in vessels over 8mm in diameter.

The BIB® PTA Dilatation Catheter is recommended for Percutaneous Transluminal Angioplasty (PTA) of the femoral, iliac and renal arteries.

These catheters are not designed to be used in the coronary arteries.

Refer to the Instructions for Use for relevant warnings, precautions, complications, and contraindications.

Rx only

BIB and CP Stent are registered trademarks of NuMED, Inc.

©2021 B. Braun Interventional Systems Inc.

CV-2091 8/21


Distributed by:

B. Braun Interventional Systems Inc.

Part of the B. Braun Group of Companies

Bethlehem, PA 18018 | USA

Tel 877-836-2228 | Fax 610-849-1334 | www.bisusa.com

 Not made with natural rubber latex

 Not made with DEHP