




Designed to Perform
at the Volumes and
Pressures You Demand

INFLATION DEVICES

Precise Control with Repeatable, Accurate
Pressure for Use in The Most Demanding
Interventional Procedures

Winged Locking Inflation Device



25cc/30ATM Winged Locking

REF Number	Product Description	Case Quantity
622510	25cc/30ATM Inflation Device with 3-way Stopcock	10
622511	25cc/30ATM Inflation Device Kit with 3-way stopcock, single Y-Connector, Metal Insertion Tool, and LTD® Torque Device	10

The 25cc/30ATM barrel and pressure gauge accommodate most procedures.

25cc clear polycarbonate syringe barrel provides rapid inflation and deflation and permits easy debubbling.


Winged locking mechanism permits easy unlocking and deflation.

Solace™ QL® Press-to-Release Inflation Device

20cc/30ATM Solace™ QL®

REF Number	Product Description	Case Quantity
622517	20cc/30ATM Inflation Device with 3-way Stopcock	10

The Solace™ QL®, Press-to-Release Inflation Device combines enhanced design features to provide a high quality, economical solution for balloon catheter dilation and fluid dispensing. The Solace™ QL® features a spring loaded locking mechanism which allows users to easily prep the device and quickly inflate or deflate a balloon, even at high pressures.



SwitchBack™ QL® Inflation Device

The SwitchBack™ QL® Inflation Devices with luminescent analog pressure gauge are available for a wide array of clinical applications such as balloon dilation and stent delivery.

Low-profile design and ergonomic handle make them easy to manipulate. The Quick-Latch mechanism allows for one-handed deflation, even when fully pressurized.

The clear blue polycarbonate syringe barrel accommodates small to very large balloon sizes and provides easy viewing of syringe contents and purging of air bubbles.



14cc/55ATM SwitchBack™ QL®

REF Number	Product Description	Case Quantity
622520	14cc/55ATM Inflation Device with 3-way Stopcock	10

25cc/40ATM SwitchBack™ QL®

REF Number	Product Description	Case Quantity
622518	25cc/40ATM Inflation Device with 3-way Stopcock	10

40cc/15ATM SwitchBack™ QL® Voluminous™

REF Number	Product Description	Case Quantity
622513	40cc/15ATM Inflation Device with 3-way Stopcock	10

60cc/15ATM SwitchBack™ QL® Voluminous™

REF Number	Product Description	Case Quantity
622519	60cc/15ATM Inflation Device with 3-way Stopcock	10

Inflation Device Accessories:

REF Number	Product Description	Case Quantity
610240	Metal Guidewire Insertion Tool	10
610400	Single (9.5F Lumen) Y-Connector Kit with Metal Insertion Tool	5
610401	Double (9.5F Lumens) Y-Connector Kit with Metal Insertion Tool	5
610420	Single (9.5F Lumen) Y-Connector Kit with Metal Insertion Tool and LTD® Torque Device	5



610240



610400



610401



610420

Indications for Use:

The Metal Guidewire Insertion Tool is designed to facilitate guidewire entry into the Y-Connector.

The Single and Double Y-Connectors provide a method for aspiration and flushing during performance of dilation.

The Torque Device is used to facilitate the manipulation of a guidewire into a vein or artery.

The 25cc/30ATM Winged Locking Inflation Device is recommended for use during the performance of balloon angioplasty procedures to inflate and deflate the balloon and to measure the pressure within the balloon during the procedure.

The Switchback™ QL® and Solace QL® Inflation Devices are recommended for use while performing balloon dilation procedures to inflate the balloon, monitor the pressure within the balloon and deflate the balloon.

Refer to the Instructions for Use for relevant warnings, precautions, complications, and contraindications. These devices are designed for single use only.

Rx only

CV-967 Rev. 9/19

LTD is a registered trademark of B. Braun Medical, Inc.

QL is a registered trademark of Atrion Medical Products Inc.

Switchback is a trademark of Atrion Medical Products Inc.

Voluminous is a trademark of B. Braun Interventional Systems Inc.

©2019 B. Braun Interventional Systems Inc.

~~PVC~~ Not made with PVC
~~LATEX~~ Not made with natural rubber latex
~~DEHP~~ Not made with DEHP

Distributed by:

B. Braun Interventional Systems Inc. | Part of the B. Braun Group of Companies

824 Twelfth Avenue | Bethlehem, PA 18018 | USA

Tel 877-836-2228 | Fax 610-849-1334 | www.bisusa.org

