

0.018in
15° Angled
and Fixed
Straight Tip

FOR YOUR ADVANCED AND
COMPLEX INTERVENTIONAL
PROCEDURES

Vascular Access Guidewires

Angled and Fixed Straight Guidewires
Available in Various Guidewire Tip Materials and Lengths

Vascular Access Guidewires

Platinum and Stainless Steel Tip Fixed Straight

Reference Number	Guidewire Diameter	Guidewire Length	Guidewire Material	Guidewire Tip	Case Qty.
610410	0.018in	40cm	NT Mandrel	PT Coil Fixed Straight	10
610411	0.018in	40cm	NT Mandrel	SS Coil Fixed Straight	10
610412	0.018in	40cm	SS Mandrel	PT Coil Fixed Straight	10
610413	0.018in	40cm	SS Mandrel	SS Coil Fixed Straight	10
610414	0.018in	80cm	NT Mandrel	PT Coil Fixed Straight	10
610415	0.018in	80cm	NT Mandrel	SS Coil Fixed Straight	10
610416	0.018in	80cm	SS Mandrel	PT Coil Fixed Straight	10
610417	0.018in	80cm	SS Mandrel	SS Coil Fixed Straight	10

Gold Tip 15° Angled and Fixed Straight

Reference Number	Guidewire Diameter	Guidewire Length	Guidewire Material	Guidewire Tip	Case Qty.
612133	0.018in	60cm	SS Mandrel	AU Coil Fixed Straight	10
612134	0.018in	40cm	NT Mandrel	AU Coil Fixed Straight	10
612135	0.018in	60cm	NT Mandrel	AU Coil Fixed Straight	10
612136	0.018in	80cm	NT Mandrel	AU Coil Fixed Straight	10
612150	0.018in	180cm	NT Mandrel	AU Coil Fixed Straight	10
612151	0.018in	80cm	NT Mandrel	AU Coil 15° Angled	10
612152	0.018in	180cm	NT Mandrel	AU Coil 15° Angled	10
612153	0.018in	300cm	NT Mandrel	AU Coil Fixed Straight	10
612154	0.018in	300cm	NT Mandrel	AU Coil 15° Angled	10

Indications For Use:

Guidewires are intended for use in percutaneous procedures to introduce and position catheters and other interventional devices within the coronary and peripheral vasculature. Refer to Instructions for Use for relevant warnings, precautions, complications, and contraindications.

LEGEND - Guidewire Material

NT = *Nitinol*
 SS = *Stainless Steel*
 PT = *Platinum*
 AU = *Gold*



Rx only

©2019 B. Braun Interventional Systems Inc. Part of the B. Braun Group of Companies.
 CV-2075 12/19

Distributed by:

B. Braun Interventional Systems Inc. | 824 Twelfth Avenue | Bethlehem, PA 18018 | USA
 Tel 877-836-2228 | Fax 610-849-1334 | www.bisusa.org

