



12F

14F

**FOR YOUR MOST  
COMPLEX CASES**

## **D'VILL™** Introducer

**Long-Length, Kink-Resistant Hemostasis Introducer**

For Introduction of Large Diameter  
Interventional Devices

# D'VILL™ Introducer

Available in 30cm, 65cm and 85cm lengths


## Hemostasis Valve

Designed to minimize blood loss




## Side Port and 3-Way Stopcock

For convenient flushing




## Braided Sheath

For optimal flexibility and kink resistance




## Radiopaque Marker Band

For easy imaging





Compatible with CP Stent® and BIB® Catheters



## Ordering Information:

Reference Number	Model Number	Sheath Size (Fr)	Usable Length (cm)
612841	DV1230	12	30
612846	DV1265	12	65
612844	DV1285	12	85
612842	DV1430	14	30
612847	DV1465	14	65
612845	DV1485	14	85

- Each D'VILL set includes (1) one Sheath and (1) one Dilator
- 0.035" guidewire compatibility

 Not made with natural rubber latex  Not made with DEHP

## Indications for Use:

Recommended for introduction of balloons, catheters, and other diagnostic and interventional devices into veins and/or arteries while maintaining hemostasis for a variety of diagnostic and therapeutic procedures.

Refer to the D'VILL® Instructions for Use for relevant warnings, precautions, complications, and contraindications.

## Rx only

D'VILL is a trademark of NuMED, Inc. CP Stent and BIB are registered trademarks of NuMED, Inc.  
©2019 B. Braun Interventional Systems Inc.  
CV-2060 Rev. 8/19

Distributed by:

**B. Braun Interventional Systems Inc.** | Part of the B. Braun Group of Companies

824 Twelfth Avenue | Bethlehem, PA 18018 | USA

Tel 877-836-2228 | Fax 610-849-1334 | [www.bisusa.org](http://www.bisusa.org)

