

Focus on the Intervention

Interventional Systems
B | BRAUN
SHARING EXPERTISE

Vascular Access Portfolio



Focus on the Intervention

B. Braun Interventional Systems Inc. offers a robust portfolio of access solutions that complement our interventional offerings. From micro-introducer kits designed to meet the needs of the smallest patients to long-length hemostasis introducers designed for the delivery of large-diameter stents—our line of access products allows you to focus on the intervention.

Micro-Introducer Kits

Tear Away Introducer Kits

Elite HV[®] Hemostasis Valve Introducer Kits


D'VILL[™] Long-Length Hemostasis Introducers

COMPASS[®] Guiding Introducer Sheath

Occlutech[®] Delivery Set III (ODS III)

ACCEL[®] Access Three-Piece Introducer System

Guidewires and Introducer Needles



Micro-Introducer Kits

Standard or Stiff Dilator with Multiple Needle and Guidewire Configurations

Micro-Introducer Kits 4 Fr x 10 cm

Reference Number	Sheath Size	Dilator	Needle	Guidewire Material	Guidewire Tip Material	Case Qty.
610373	4 Fr	Standard	4 cm	NT	PT	5
610374	4 Fr	Standard	4 cm	NT	SS	5
610375	4 Fr	Stiff	4 cm	NT	PT	5
613317	4 Fr	Standard	4 cm	NT	AU	5
613318	4 Fr	Stiff	4 cm	NT	AU	5
610361	4 Fr	Standard	7 cm	NT	PT	10
610362	4 Fr	Standard	7 cm	NT	SS	10
610363	4 Fr	Standard	7 cm	SS	PT	10
610364	4 Fr	Standard	7 cm	SS	SS	10
610365	4 Fr	Stiff	7 cm	NT	PT	10
610366	4 Fr	Stiff	7 cm	NT	SS	10
610390	4 Fr	Stiff	7 cm	SS	SS	10
610396	4 Fr	Stiff	7 cm	NT	PT*	10
613311	4 Fr	Standard	7 cm	NT	AU	10
613312	4 Fr	Standard	7 cm	SS	AU	10
613313	4 Fr	Stiff	7 cm	NT	AU	10

*Guidewire tip is angled

Micro-Introducer Kits 5 Fr x 10 cm

Reference Number	Sheath Size	Dilator	Needle	Guidewire Material	Guidewire Tip Material	Case Qty.
610377	5 Fr	Standard	4 cm	NT	PT	5
610378	5 Fr	Standard	4 cm	NT	SS	5
610379	5 Fr	Stiff	4 cm	NT	PT	5
610394	5 Fr	Standard	4 cm	SS	SS	5
613319	5 Fr	Standard	4 cm	NT	AU	5
613320	5 Fr	Stiff	4 cm	NT	AU	5
610367	5 Fr	Standard	7 cm	NT	PT	10
610368	5 Fr	Standard	7 cm	NT	SS	10
610369	5 Fr	Standard	7 cm	SS	PT	10
610370	5 Fr	Standard	7 cm	SS	SS	10
610371	5 Fr	Stiff	7 cm	NT	PT	10
610372	5 Fr	Stiff	7 cm	NT	SS	10
610391	5 Fr	Stiff	7 cm	SS	SS	10
613314	5 Fr	Standard	7 cm	NT	AU	10
613315	5 Fr	Standard	7 cm	SS	AU	10
613316	5 Fr	Stiff	7 cm	NT	AU	10

Seamless™ Micro-Introducer Kits 4 Fr x 5 cm

Reference Number	Sheath Size	Dilator	Needle Non-echogenic	Guidewire Material	Guidewire Tip Material	Case Qty.
614713	4 Fr x 5 cm	Standard	2.5 cm	NT	AU	5
614714	4 Fr x 5 cm	Stiff	2.5 cm	NT	AU	5

Micro-Introducer Kit Contents:

One micro-introducer
One 21-gauge echogenic needle
One 0.018 in. x 40 cm guidewire



Stiffen Dilator

Reference Number	French Size	Length	Dilator	Case Qty.
610392	3 Fr	12.5 cm	Stiff	5

Guidewire Material:

SS = Stainless Steel NT = Nitinol AU = Gold PT = Platinum


Tear Away Introducer Kits

Tear Away Introducer Kits

Reference Number	Tear Away Introducer	Case Qty.
612601	5 Fr x 14 cm	5
612602	6 Fr x 14 cm	5
612603	7 Fr x 14 cm	5
612604	8 Fr x 14 cm	5
612606	9 Fr x 14 cm	5
612608	10 Fr x 14 cm	5
612609	10.5 Fr x 14 cm	5
612610	11 Fr x 14 cm	5
612611	12 Fr x 14 cm	5
612612	13 Fr x 14 cm	5
612613	14 Fr x 14 cm	5
612614	15 Fr x 14 cm	5
612615	16 Fr x 14 cm	5

Tear Away Introducer Kit Contents:

One tear away introducer, PTFE
 One 18-gauge x 7 cm introducer XTW (extra-thin wall) needle
 One 10 mL syringe
 One 0.038 in. x 50 cm double flex tip/3 mm J-straight, stainless steel guidewire



~~LATEX~~
 Not made with natural rubber latex

~~DEHP~~
 Not made with DEHP

~~PVC~~
 Not made with PVC

Mini Tear Away Introducer Kits

Reference Number	Mini Tear Away Introducer	Case Qty.
612645	4 Fr x 7 cm	10
612646	5 Fr x 7 cm	10
612647	6 Fr x 7 cm	10
612648	4 Fr x 10 cm	10
612649	5 Fr x 10 cm	10
612650	6 Fr x 10 cm	10

Mini Tear Away Introducer Kit Contents:

One tear away introducer, PTFE
 One 21 gauge x 7 cm introducer XTW needle
 One 0.018 in. x 45 cm Nitinol Mandrel, platinum tip guidewire

Stand-Alone Tear Away Introducer

Reference Number	Tear Away Introducer	Case Qty.
612780*	4 Fr x 14 cm	10
612781	5 Fr x 14 cm	10
612782	6 Fr x 14 cm	10
612783	7 Fr x 14 cm	10
612616	8 Fr x 14 cm	10
612617	9 Fr x 14 cm	10
612619	10 Fr x 14 cm	10
612621	11 Fr x 14 cm	10
612622	12 Fr x 14 cm	10
612624	14 Fr x 14 cm	10
612626	16 Fr x 14 cm	10

Recommended guidewire: 0.038 in.

*Recommended guidewire: 0.028 in.

Stand-Alone Tear Away Long Introducer

Reference Number	Tear Away Introducer	Case Qty.
612791	5 Fr x 30 cm	5
612792	6 Fr x 30 cm	5
612793	7 Fr x 30 cm	5
612627	8 Fr x 30 cm	5
612628	9 Fr x 30 cm	5
612629	10 Fr x 30 cm	5
612630	11 Fr x 30 cm	5
612631	12 Fr x 30 cm	5
612632	13 Fr x 30 cm	5
612633	14 Fr x 30 cm	5
612634	15 Fr x 30 cm	5
612635	16 Fr x 30 cm	5
612852	6 Fr x 44 cm	5
612637	9 Fr x 44 cm	5
612640	12 Fr x 44 cm	5

Recommended guidewire: 0.038 in.



Elite HV[®]

Hemostasis Valve Introducer Kits

Reference Number	French Size	Introducer Sheath Length	Guidewire Shaft Material	Guidewire Tip Material	Guidewire Diameter and Length	Case Qty.
612930	3 Fr	7 cm	NT	PT/SS*	0.018 in x 40 cm	10
612800	3 Fr	7 cm	NT	AU/SS*	0.018 in x 40 cm	10
612801	3 Fr	7 cm	SS	SS	0.021 in x 45 cm	10
612802	4 Fr	7 cm	SS	SS	0.021 in x 45 cm	10
612803	5 Fr	7 cm	SS	SS	0.021 in x 45 cm	10
612804	6 Fr	7 cm	SS	SS	0.021 in x 45 cm	10
612805	7 Fr	7 cm	SS	SS	0.021 in x 45 cm	10

*This kit contains a dual floppy-tipped guidewire

Guidewire Material:

SS = Stainless Steel NT = Nitinol AU = Gold PT = Platinum



Hemostasis Valve Introducer Kit Contents:

One Elite HV Introducer Sheath
One guidewire



Not made with
natural rubber latex



Not made
with DEHP

D'VILL[™] Introducers

Long-Length, Kink-Resistant Hemostasis Introducer

Reference Number	Sheath Size	Usable Length
612841	12 Fr	30 cm
612846	12 Fr	65 cm
612844	12 Fr	85 cm
612842	14 Fr	30 cm
612847	14 Fr	65 cm

0.035 in. guidewire compatibility




Not made with
natural rubber latex



Not made
with DEHP

COMPASS® Guiding Introducer Sheath

Kink-Resistant Introducer Sheath



Size	Length		Straight	Multipurpose MP	Left Internal Mammary Artery LIMA	Renal Double Curve RDC	Hockey Stick HS	Inner Diameter (minimum)	Dilator Guidewire Compatibility	Case Qty.
4 Fr	15 cm*	Red	613944**					1.60 mm / 0.063 in	0.035 in and 0.018 in	5
4 Fr	45 cm		613945					1.60 mm / 0.063 in	0.035 in	5
4 Fr	60 cm		613946**					1.60 mm / 0.063 in	0.035 in	5
4 Fr	90 cm		613947					1.60 mm / 0.063 in	0.035 in	5
5 Fr	15 cm*	Grey	613948**					1.92 mm / 0.075 in	0.035 in and 0.018 in	5
5 Fr	45 cm		613949	613950	613972	613973	613974	1.92 mm / 0.075 in	0.035 in	5
5 Fr	60 cm		613951**					1.92 mm / 0.075 in	0.035 in	5
5 Fr	90 cm		613952**	613953				1.92 mm / 0.075 in	0.035 in	5
6 Fr	15 cm*	Green	613954**					2.25 mm / 0.088 in	0.035 in and 0.018 in	5
6 Fr	45 cm		613955	613956	613975**	613976	613977**	2.25 mm / 0.088 in	0.035 in	5
6 Fr	60 cm		613957					2.25 mm / 0.088 in	0.035 in	5
6 Fr	90 cm		613958**	613959**				2.25 mm / 0.088 in	0.035 in	5
7 Fr	15 cm*	Orange	613960					2.55 mm / 0.100 in	0.035 in and 0.018 in	5
7 Fr	45 cm		613961	613962	613978**	613979**	613980**	2.55 mm / 0.100 in	0.038 in	5
7 Fr	60 cm		613963**					2.55 mm / 0.100 in	0.038 in	5
7 Fr	90 cm		613964	613965				2.55 mm / 0.100 in	0.038 in	5
8 Fr	15 cm*	Blue	613966**					2.88 mm / 0.113 in	0.035 in and 0.018 in	5
8 Fr	45 cm		613967	613968**				2.88 mm / 0.113 in	0.038 in	5
8 Fr	60 cm		613969**					2.88 mm / 0.113 in	0.038 in	5
8 Fr	90 cm		613970**	613971**				2.88 mm / 0.113 in	0.038 in	5


*15 cm versions come with two dilators (0.035 in. and 0.018 in.)

** Made to order



Occlutech® Delivery Set III

ODS III



Set Contents:

- (A) One dilator
- (B) One sheath
- (C) One loader

Sheath Type	Reference Number	Size	Usable Length Sheath	Usable Length Loader	Distal Angle	French Size Color Code
45° Regular Length 	98DS006	6 Fr	79 cm	8.5 cm	45°	
	98DS007	7 Fr	79 cm	8.5 cm	45°	
	98DS008	8 Fr	79 cm	12.5 cm	45°	
	98DS009	9 Fr	79 cm	12.5 cm	45°	
	98DS010	10 Fr	79 cm	12.5 cm	45°	
	98DS011	11 Fr	79 cm	16.5 cm	45°	
	98DS012	12 Fr	79 cm	16.5 cm	45°	
	98DS014*	14 Fr	79 cm	16.5 cm	45°	
45° Long Length 	98DL006*	6 Fr	110 cm	8.5 cm	45°	
	98DL007*	7 Fr	110 cm	8.5 cm	45°	
	98DL008*	8 Fr	110 cm	12.5 cm	45°	
	98DL009*	9 Fr	110 cm	12.5 cm	45°	
	98DL010*	10 Fr	110 cm	12.5 cm	45°	
180° Regular Length 	98US006*	6 Fr	79 cm	8.5 cm	180°	
	98US007*	7 Fr	79 cm	8.5 cm	180°	
	98US008*	8 Fr	79 cm	12.5 cm	180°	
	98US009*	9 Fr	79 cm	12.5 cm	180°	
	98US010*	10 Fr	79 cm	12.5 cm	180°	
	98US011*	11 Fr	79 cm	16.5 cm	180°	

 Not made with natural rubber latex

 Not made with DEHP

 Not made with PVC

* Made to order

ACCEL® Access

Three-Piece Introducer System with TrueGlide® Hydrophilic Coating

ACCEL Access Three-Piece Introducer

Reference Number	Introducer Size	Introducer Length	Case Qty.
612137	6 Fr	20 cm	10

Single Echogenic Trocar or Chiba Introducer Needle with Depth Markings

Reference Number	Needle Type	Needle Gauge	Needle Length	Case Qty.
612130	Trocar	21 Ga	15 cm	5
612131	Chiba	22 Ga	15 cm	5
612132	Chiba	22 Ga	20 cm	5

ACCEL Access Three-Piece Introducer with Echogenic Trocar or Chiba Needle Set

Reference Number	Introducer	Needle Type	Needle Gauge	Needle Length	Case Qty.
612144	6 Fr x 20 cm	Chiba	22 Ga	15 cm	10
612145	6 Fr x 20 cm	Trocar	21 Ga	15 cm	10

ACCEL Access Three-Piece Introducer with Echogenic Trocar Needle and Guidewire Set

Reference Number	Introducer	Needle Type	Needle Gauge	Needle Length	Guidewire Diameter	Guidewire Length	Guidewire Material	Guidewire Tip	Case Qty.
612138	6 Fr x 20 cm	Trocar	21 Ga	15 cm	0.018 in	60 cm	NT Mandrel	AU Coil	10
612139	6 Fr x 20 cm	Trocar	21 Ga	15 cm	0.018 in 0.038 in	60 cm 145 cm	SS Mandrel PTFE	AU Coil 3 mm J-Straight	5
612140	6 Fr x 20 cm	Trocar	21 Ga	15 cm	0.018 in 0.038 in	60 cm 145 cm	NT Mandrel PTFE	AU Coil 3 mm J-Straight	5
612141	6 Fr x 20 cm	Trocar	21 Ga	15 cm	0.018 in	60 cm	SS Mandrel	AU Coil	10

ACCEL Access Three-Piece Introducer with Echogenic Chiba Needle and Guidewire Set

Reference Number	Introducer	Needle Type	Needle Gauge	Needle Length	Guidewire Diameter	Guidewire Length	Guidewire Material	Guidewire Tip	Case Qty.
612142	6 Fr x 20 cm	Chiba	22 Ga	15 cm	0.018 in	60 cm	NT Mandrel	AU Coil	10
612143	6 Fr x 20 cm	Chiba	22 Ga	15 cm	0.018 in	60 cm	SS Mandrel	AU Coil	10

Guidewire Material:

SS = Stainless Steel NT = Nitinol AU = Gold



Not made with
natural rubber latex




Not made
with DEHP



Not made
with PVC



Vascular Access Guidewires




Platinum and Stainless Steel Fixed Straight Tip

Reference Number	Guidewire Diameter	Guidewire Length	Guidewire Material	Flexible Guidewire Tip, 5 cm	Case Qty.
610410	0.018 in	40 cm	NT Mandrel	PT Coil	10
610411	0.018 in	40 cm	NT Mandrel	SS Coil	10
610412	0.018 in	40 cm	SS Mandrel	PT Coil	10
610413	0.018 in	40 cm	SS Mandrel	SS Coil	10
610414	0.018 in	80 cm	NT Mandrel	PT Coil	10
610415	0.018 in	80 cm	NT Mandrel	SS Coil	10
610416	0.018 in	80 cm	SS Mandrel	PT Coil	10
610417	0.018 in	80 cm	SS Mandrel	SS Coil	10

Gold Tip Fixed Straight and 15° Angled Tip

Reference Number	Guidewire Diameter	Guidewire Length	Guidewire Material	Flexible Guidewire Tip, 5 cm	Case Qty.
612133	0.018 in	60 cm	SS Mandrel	AU Coil Fixed Straight	10
612134	0.018 in	40 cm	NT Mandrel	AU Coil Fixed Straight	10
612135	0.018 in	60 cm	NT Mandrel	AU Coil Fixed Straight	10
612136	0.018 in	80 cm	NT Mandrel	AU Coil Fixed Straight	10
612150	0.018 in	180 cm	NT Mandrel	AU Coil Fixed Straight	10
612151	0.018 in	80 cm	NT Mandrel	AU Coil 15° Angled	10
612152	0.018 in	180 cm	NT Mandrel	AU Coil 15° Angled	10
612153	0.018 in	300 cm	NT Mandrel	AU Coil Fixed Straight	10
612154	0.018 in	300 cm	NT Mandrel	AU Coil 15° Angled	10



Guidewire Material:

SS = Stainless Steel NT = Nitinol AU = Gold PT = Platinum


~~LATEX~~
Not made with
natural rubber latex

~~DEHP~~
Not made
with DEHP

~~PVC~~
Not made
with PVC

~~BPA~~
Not made with
bisphenol A

Introducer Needles



Safecan™ Safety Introducer Needles

Reference Number	Size	Needle Length	Hub Color	Needle Bevel	Echogenic Needle Tip	Guidewire Compatibility	Case Qty.
613444	21 Ga	7 cm	Clear	A	No	0.018 in	25
613445	21 Ga	7 cm	Clear	A	Yes	0.018 in	25



Not made with
natural rubber latex



Not made
with DEHP



Not made
with PVC

Extra Sharp Introducer Needles

Reference Number	21 Gauge Needle	Hub Color	Needle Bevel	Guidewire Compatibility	Case Qty.
614820	4 cm Echogenic	Green	B	≤ 0.021 in	25
614821	7 cm Echogenic	Green	B	≤ 0.021 in	25



Not made with
natural rubber latex



Not made
with DEHP



Not made
with PVC

Standard Introducer Needles

Reference Number	21 Gauge Needle	Hub Color	Needle Bevel	Guidewire Compatibility	Case Qty.
614702	2.5 cm Non-echogenic	Green	B	≤ 0.021 in	25
613100	4 cm Echogenic	Green	B	≤ 0.021 in	25
613101	7 cm Echogenic	Green	B	≤ 0.021 in	25
613102	15 cm Echogenic	Green	B	≤ 0.021 in	25




Not made with
natural rubber latex



Not made
with DEHP



Not made
with PVC




Vascular Access Products

Introducers, Guidewires and Needles




Micro-Introducer Kits

- Micro-introducer sheaths and dilators available in 4 Fr and 5 Fr sizes, 5 cm and 10 cm lengths, standard and stiff
- Seamless dilator-to-sheath transition facilitates smooth passage into the vessel
- 2.5 cm, 4 cm and 7 cm needles and several guidewire configurations available for varying clinical needs



Tear Away Introducer Kits

- Available in 4-16 Fr sizes in 7 cm, 10 cm, 14 cm, 30 cm and 44 cm lengths
- Kink-resistant PTFE sheath material assures sheath integrity and eases device insertion
- Seamless dilator-to-sheath transition facilitates atraumatic device insertion




Hemostasis Valve Introducer Kits

- Available in 3-7 Fr sizes in 7 cm lengths
- Designed for easy catheter movement through the hemostasis valve during multiple catheter insertions
- Seamless dilator-to-sheath transition for atraumatic access



D'VILL™ Introducers

- For introduction of large diameter interventional devices, available in 12 Fr and 14 Fr sizes in 30 cm, 65 cm and 85 cm lengths
- Braided sheath designed to provide optimal flexibility and kink resistance



COMPASS®

- Available in 4-8 Fr sizes and 15 cm, 45 cm, 60 cm and 90 cm lengths
- Straight, multipurpose, hockey stick, RDC and LIMA tip options



Occlutech® Delivery Set III (ODS III)

- Available in 6-14 Fr sizes, 79 cm and 110 cm lengths
- 45- and 180-degree-shaped tip configurations



ACCEL® Access Three-Piece Introducers

- 6 Fr Introducer with TrueGlide™ Hydrophilic Sheath Coating for smooth insertion
- For vascular and non-vascular applications




Vascular Access Guidewires

- 0.018 in. Fixed straight and angled tip guidewires, available in various lengths and guidewire tip materials



Safecan™ Safety Introducer Needles

- 18 Ga. and 21 Ga. Safety Introducer Needles in 7 cm length, echogenic and non-echogenic tips
- Designed to provide clinician and patient safety and reduce the risk of needlestick injury



Extra Sharp and Standard Introducer Needles

- 21 Ga. Extra Sharp Echogenic Introducer Needles in 4 cm and 7 cm lengths
- 21 Ga. Standard Introducer Needles in 2.5 cm, 4 cm, 7 cm and 15 cm lengths, echogenic and non-echogenic tips

Indications For Use:

Micro Introducer Kits are intended to introduce up to a 0.038 in. guidewire or catheter into the vascular system following a small-gauge needle stick.

Coaxial Dilator Sets are intended to introduce up to a 0.038 in. guidewire or catheter into the vascular system following a small-gauge needle stick.

Tear Away Introducers are used for the percutaneous introduction of diagnostic or therapeutic devices, such as catheters and pacing leads into the vasculature.

Elite HV® Hemostasis Valve Introducer System is indicated for use in percutaneous procedures to introduce catheters and other intravascular devices into the vasculature.

D'VILL™ Introducers are recommended for introduction of balloons, catheters and other diagnostic and interventional devices into veins and/or arteries while maintaining hemostasis for a variety of diagnostic and therapeutic procedures.

Compass® Guiding Introducer Sheath is indicated to be used for introduction of interventional and diagnostic devices into the peripheral (and coronary) vasculature. The device is also intended to be used within a pediatric population.

Occlutech® Delivery Set III is intended for the introduction of diagnostic and therapeutic devices into the human vasculature, including but not limited to intracardiac, renal or other placements.

ACCEL® Three-Piece Introducer System is indicated for use in percutaneous introduction and placement of catheters and guidewires.

Trocar and Chiba Introducer needles are used for percutaneous introduction of guidewires.

Guidewires are intended for use in percutaneous procedures to introduce and position catheters and other interventional devices within the coronary and peripheral vasculature.

Needles are used for the percutaneous introduction of guidewires.

Refer to the Instructions For Use for relevant warnings, precautions, complications and contraindications.

Elite HV is a registered trademark of Galt Medical Corp.

D'VILL is a trademark of NuMED, Inc.

Safecan is a trademark of B. Braun Medical Inc.

ACCEL, TrueGlide and COMPASS are registered trademarks of B. Braun Interventional Systems Inc.

Seamless is a trademark of B. Braun Interventional Systems Inc.

Occlutech is a registered trademark of Occlutech Holding AG.

For more information

Rx only.

©2026 B. Braun Interventional Systems Inc.

Distributed by:

B. Braun Interventional Systems Inc. | a member of the B. Braun Group of Companies

Bethlehem, PA 18018 | USA | 877-836-2228 | bisusa.com

